

## Micro B3 EU-US Bioinformatics Course: Day 2

### Command line BLAST, parsing output

#### 1. Go back to the home directory and start an interactive torque job

```
$ cd ~  
$ int
```

#### 2. Create a few output directories

```
$ mkdir output_repository  
$ cd output_repository  
$ mkdir blastn  
$ mkdir blastp  
$ mkdir tblastn  
$ mkdir extracted  
$ cd ~
```

#### 3. Performing a BLASTn

```
$ blastn -db ~/microb3-shared/blast/databases/small.nuc  
-query ~/microb3-shared/blast/queries/small.nuc.read.fasta  
-out ~/output_repository/blastn/small.raw -outfmt 0
```

#### 4. Checking out the results

```
$ cd output_repository/blastn  
$ less small.raw
```

#### 5. Another format for output (more compact format)

```
$ cd ~  
$ blastn -db ~/microb3-shared/blast/databases/small.nuc  
-query ~/microb3-shared/blast/queries/small.nuc.read.fasta  
-out ~/output_repository/blastn/small.btab -outfmt 6
```

#### 6. Checking out the results

```
$ cd output_repository/blastn  
$ less small.btab
```

#### 7. Performing a BLASTp

```
$ cd ~
$ blastp -query ~/microb3-shared/blast/queries/small.pep.orf
.fasta -db ~/microb3-shared/blast/databases/small.pep
-out output_repository/blastp/small.raw -outfmt 0
```

## 8. Checking out the results

```
$ cd output_repository/blastp
$ less small.raw
```

## 9. Bigger BLASTp

```
$ cd ~
$ blastp -query ~/microb3-shared/blast/queries/dna_polymerase_I.fasta
-db ~/microb3-shared/blast/databases/cbay.orf.pep
-out ~/output_repository/blastp/dnapola.btab -outfmt 6
```

## 10. Checking out the results

```
$ cd output_repository/blastp
$ less dnapola.btab
```

## 11. How many hits?

```
$ wc -l dnapola.btab
```

## 12. How many unique query sequences, unique subject sequences?

```
$ cut -f1 dnapola.btab | sort -u | wc -l
```

```
$ cut -f2 dnapola.btab | sort -u | wc -l
```

## 13. Extracting query sequences

```
$ extract_blast_hits.pl -m
```

```
$ extract_blast_hits.pl --btab ~/output_repository/blastp/dnapola.btab
--fasta ~/microb3-shared/blast/databases/cbay.orf.pep
--out ~/output_repository/extract.subj.fasta --subj
```

## 14. Using other BLAST flavors (tBLASTn)

```
$ cd ~
```

```
$ tblastn -db ~/microb3-shared/blast/databases/cbay.reads -  
query ~/microb3-shared/blast/queries/dna_polymerase_I.fasta -  
out ~/output_repository/tblastn/cbay.dnapol.btab -outfmt 6
```

```
$ tblastn -db ~/microb3-shared/blast/databases/cbay.reads -  
query ~/microb3-shared/blast/queries/dna_polymerase_I.fasta -  
out ~/output_repository/tblastn/cbay.dnapol.raw -outfmt 0
```

## 15. Tinkering with word size

```
$ exit
```

```
$ cd ~
```

```
$ less ~/microb3-shared/blast/queries/crispr.spacers.gos.fasta
```

```
$ less ~/microb3-shared/bin/blastn-crispr.qs
```

```
$ qsub ~/microb3-shared/bin/blastn-crispr.qs
```

```
$ wc -l ~/output_repository/blastn/crispr.btab
```

```
$ less ~/microb3-shared/bin/blastn-crispr-word.qs
```

```
$ qsub ~/microb3-shared/bin/blastn-crispr-word.qs
```

```
$ wc -l ~/output_repository/blastn/crispr.7.btab
```